NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • http://www.nti.de

DE Gebrauchsanweisung RPR System

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RPR Root post and core system

Product description:

The RPR system is a passive root post and core system with conical root posts. The post material is pure titanium.

Plastic posts are also available in the RPR root post and core system for customised post and cores.

Indication/Intended purpose:

For building up broken teeth or teeth severely destroyed by caries or trauma.

Contraindication:

Teeth with damage that extends below the gingival margin. Broken roots or roots with a fracture or crack.

Information on cleaning, disinfection and sterilisation of titanium root posts and accessories:

The titanium root posts and instruments of the system can be disinfected and sterilised in an autoclave.

Please note: Use only approved disinfectants without alcohol. After disinfection, rinse the root posts and instruments under running water, dry and shrink-wrap in foil to protect against recontamination. Sterilisation is completed at 134°C, for 5 min.

(Detailed preparation recommendations according to DIN EN ISO 17664 ☞ www.nti.de).

Information on preparing plastic posts:

The plastic posts cannot be disinfected or sterilised.

Medical requirements for the use of prefabricated RPR root posts:

1. The root must be periodontally intact and caries-free. Root treatment must be successfully completed.
2. The lesion should not extend below the gingival margin.
3. A ferrule crown preparation is therefore possible with a minimum of 1.5 - 2.0 mm in the tooth structure.

All instruments and posts of the RPR system should be secured to prevent aspiration and swallowing. The use of a rubber dam is recommended.

Only use conditioned RPR titanium root posts and instruments of the RPR system.

Use of the titanium root posts:

Please note: Titanium root posts are intended for single use. Do not reuse!

1. Smooth the surface of the tooth using a diamond instrument.
2. Select as large a root post as possible on the X-ray with the aid of the template. The root post should finish 2 mm before the apical end of the foramen. It must fit flush to the walls of the canal for a minimum depth of 2/3 of the canal length.
3. If it is not possible to complete point 2, a customised post and core should be fabricated.
4. The wall thickness of the root should not be less than 1.5 mm, to ensure sufficient stability.
5. The post should not be shortened.
6. The root posts have a head. A root facer produces a perfect rest seat for the post head on the tooth. Drill the rest seat in the root surface using a root facer, which is the same size as the post. Motor speed: opt. 1,000-1,500 rpm. Do not tilt the cutter (risk of fracture)!
7. Attach the depth stop, which corresponds to the length of the post, onto the predrill.
8. Extend the root canal using a handpiece with the appropriate predrill.
Motor speed: opt. 500-1,000 rpm.

9. Provide the root canal with retention using the roughening instrument. Secure the instrument against aspiration and swallowing. Roughening is not necessary, if an adhesive system is to be used.
10. Clean the root canal with water and dry with paper points.
11. The posts can be luted using phosphate cement. Adhesive systems are also suitable. Always follow the respective manufacturer's instructions.
12. Coat the posts with the selected luting system and insert into the root canal. Secure against aspiration and swallowing, for example using dental floss.
13. The post can also be vibrated into the root canal using an ultrasonic probe.
14. Remove the cement residue after curing.
15. Fabricate the core using composite or cement.
16. The crown preparation should extend a minimum of 1.5-2.0 mm into the healthy tooth structure. This is the only way to ensure that the forces produced are transferred to the residual tooth.

Use of the plastic posts:

Please note: Plastic posts were developed for single use only. Do not reuse!

There are two options for using the plastic posts

1. The same medical requirements apply as for titanium posts.
2. Smooth the surface of the tooth using a diamond instrument.
3. Select as large a root post as possible on the X-ray with the aid of the template. The root post should finish 2 mm before the apical end of the foramen. It must fit flush to the walls of the canal for a minimum depth of 2/3 of the canal length.
4. If it is not possible to complete point 3, a customised post and core should be fabricated.
5. The wall thickness of the root should not be less than 1.5 mm, to ensure sufficient stability.
6. The post can be shortened incisally/occlusally.
7. Extend the root canal using a handpiece with the appropriate predrill. Motor speed: opt. 500-1,000 rpm

The procedure is now different according to application:

Indirect core – recommended technique:

1. Insert the selected post in the canal and take an impression over it.
2. Send the impression with the post to the laboratory.
3. The laboratory fabricates a cast post and core.

Direct core:

1. Separate the root canal using Vaseline or silicone.
2. Insert the post into the canal and fabricate the core using a resin that is approved for intraoral use.
3. Prepare the core.
4. Send the post and core to the laboratory.
5. The laboratory casts the post and core.

In both techniques, cement the post and core, complete the preparation and take an impression.

General information:

The RPR root post systems are developed and approved only for use in the dental sector. The system should be used according to the instructions for use.

The instructions for use were compiled to the best of our knowledge and in good faith and according to scientific methods at the time of publishing.

Any other use of the product other than that given in the product information is the sole responsibility of the user.

The RPR root post system is approved only for use by dentists.

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