
NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • www.nti.de

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GB Instructions for use for conical FiberMaster glass fibre root posts

The conical FiberMaster glass fibre posts may only be used on patients by dentists.

The conical FiberMaster glass fibre posts have been developed to stabilise and restore teeth with partially damaged crowns using a composite. Please read these instructions carefully prior to using the conical FiberMaster posts.

The FiberMaster material has been specially designed for dental applications and must be used according to the recommendations for use.

The manufacturer does not accept any liability for damage caused by failure to observe these instructions or incorrect use. The user is personally responsible for checking the suitability of the material for the intended purpose. This is particularly necessary if the intended application is not included in the instructions for use.

Contents:

1. Conical FiberMaster indications

1.1. Prerequisites for the use of FiberMaster posts

2. Recommendations for use

3. Information on preparation

1. Conical FiberMaster indications

Teeth with badly damaged crowns and fractured teeth with no other root lesions and which are still completely periodontally intact. Teeth with fractured or cracked roots are contraindicated.

1.1. Prerequisites for the use of the FiberMaster system:

- Root treatment must have been successfully completed.
- The root must be periodontally intact and caries-free.
- The lesion should not extend below the gingival margin.
- There must be continuous contact between the wall and the post along at least 2/3 of the canal length.
- Sufficient root wall thickness (at least 1.5 mm).
- As large a post diameter as possible.
- The occlusal or incisal surface of the stump must be prepared for an abutment.
- A ferrule preparation at least 1.5 to 2.0 mm into the natural hard substance is required. The posts may fracture without a ferrule preparation. In such cases, the creation of customised metal abutments is recommended.
- All instruments and FiberMaster system posts are to be secured to prevent aspiration and swallowing by the patient. The use of a rubber dam is recommended.
- All FiberMaster system instruments must only be used once they have been prepared. The posts must be disinfected with 70% medical alcohol and all traces of grease must have been cleaned off (see point 3).

2. Recommendations for use:

- The FiberMaster posts are designed for single use. Do not re-use! Smooth the surface of the tooth with a diamond instrument.
- Select the size of the root post using the x-ray. It must be decided at this stage whether a post with or without a head is to be used. In the case of teeth with several roots, multiple posts with or without a head can be used depending on the space available.
- Remove the root filling.
- With the FiberMaster posts with a head, a root facer

is used to prepare the stump surface for positioning of the head and cut to a depth of approx. 1 mm (speed opt. 1,000 - 1,500 rpm). This stage is not required for the conical FiberMaster posts without a head.

- Fit the correct depth gauge to the drill.
- Enlarge the canal to the necessary depth using the drill. This is determined on the x-ray image. In order to mark the required depth, stoppers from the endodontics range can be fitted (speed: 500 - 1,000 rpm). When doing this, it must be ensured that at least 4 mm of filled canal remains in the apical area.
- Remove the conical FiberMaster from the packaging using forceps. Do not touch the posts with your fingers as oil from the skin prevents proper cementation.
- The length of the conical FiberMaster posts can be easily adapted to the prevailing conditions using a fine-grain diamond grinding disc. To this end, only use discs without holes and slits. Shorten the post extraorally. The posts are radiopaque. This means the positioning of the post can be checked on the x-ray image.
- If the posts were touched with bare hands, remove any soiling using 70% medical alcohol. Do not touch the posts with bare hands after cleaning.

Sterilisation is not possible as the temperatures would damage the material.

- An adhesive luting technique with ETCH AND RINSE must be used. A rubber dam must always be used with this technique.
- *Self-etching one-bottle systems are not suitable for cementing the glass fibre posts.
- Etch the root canal and the stump surface with 37% phosphoric acid.
- Rinse off the phosphoric acid.
- Dry using either paper tips or gently with oil-free air.
- Treat the root canal and the surface of the tooth with primer and bonder according to the manufacturer's instructions and light cure as necessary.
- If so desired, the root post can be treated with a silane. This, however, is not absolutely necessary.
- The luting composite is applied to the root post.
- A low-viscosity dual composite is to be used for cementing. High-viscosity dual composites should only be used with an ultrasonic unit. When cementing the posts, apply light pressure occlusally to the posts using a ball instrument.
- Only dual composites are to be used which also cure without light. The instructions for use issued by the manufacturer of the adhesive must always be followed! Adhesive systems should be used which are suitable for building up the abutment. The recommended systems are Luxabond and Luxacore from DMG and Panavia F2 from Kuraray.
- As with all adhesive systems, it is standard practice to use a rubber dam.

Please note:

- Do not clean the root canal with NaOCl (sodium hypochlorite) or H₂O₂ (hydrogen peroxide) as these substances can have a detrimental effect on the polymerisation of the composite. These chemicals may prevent curing and adhesion in the canal.
- Apply the adhesive to the root post and insert into the canal applying gentle pressure to allow any excess to escape.
- Use the excess to shape the abutment. The result is a monoblock made from one material. This ensures enhanced stability. Then light cure.
- **Important:** The root post does not conduct light sufficiently. That is why self-curing or dual composite adhesives are necessary.
- Make the abutment from the same material which was used to cement the post (please follow the manufacturer's instructions). The resulting monoblock offers greater stability than a block made from two different materials.
- Make a ferrule preparation at least 1.5 - 2 mm below the root surface.
- Suitable measures must be taken to secure the posts so as to prevent aspiration or swallowing by the patient during placement. Use of a rubber dam is recommended for this reason.
- If necessary, adjust the size and shape of the head intraorally using a diamond instrument. For this, a rubber dam must be used to prevent the aspiration and swallowing of glass fibres.

3. Information on preparation:

The rotary instruments are disinfected and sterilised (134°C, 5 minutes) in the normal manner before being re-used.

Detailed preparation recommendations as per DIN EN ISO 17664 www.nti.de

- The posts are disinfected with 70% medical alcohol and all traces of grease are cleaned off. Sterilisation and disinfection in a thermal disinfectant is not possible.

- Do not touch the posts with bare hands after cleaning. *Literature review Quintessenz der Zahnmedizin 7/2012 "Sind selbstkonditionierende Bondingsysteme praxistauglich? Und wenn ja wofür? Teil 2" Prof. Dr. Bernd Haller, Dr. Florian Janke. Quintessenz International July/August 2012 "Fiber post adhesion to radicular dentin: The use of acid etching prior to a one-step self-etching adhesive".

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© NTI-Kahla GmbH Rotary Dental Instruments