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Contents:

1. Product description
2. Indication for parapulpal retention pins
3. Contraindication for parapulpal retention pins
4. Notes on processing
5. Recommendations for use
- 5.1. Retopin Gold, Reto-Minipin Gold and Retopin Titanium, Reto-Minipin Titanium
- 5.2. TRI-STAR
- 5.3. TRI-JET

1. Product description

Dear customer,

Please read the instructions for use carefully before using the parapulpal retention pins and the corresponding instruments.

The parapulpal retention pins were developed for use by trained dentists and may only be used by these. The parapulpal retention pins were developed for tooth preservation. The manufacturer shall not be liable for any damages due to non-compliance with these provisions in the instructions for use or which are due to system-incompatible use.

The user is obliged to check the material for its suitability for the intended use at his/her own responsibility. This is particularly necessary if this concerns a use not mentioned in the product information.

Retopin Gold and Reto-Minipin Gold:

Pin: stainless steel (1.4305) electrolytically gold-plated

Mandrel: anodised aluminium

Retopin Titanium and Reto-Minipin Titanium:

Pin: titanium TiAl6V4

Mandrel: anodised aluminium

TRI-STAR and TRI-JET

Pin and mandrel: Stainless steel (1.4305)

2. Indication for parapulpal retention pins:

The retention pins are screwed into the dentine next to (para) the pulp.

- Additional retention for all filling materials.
- As protection against twisting for metallic root post systems.
- As retention and protection against twisting for abutments made of acrylic materials.
- For building up broken teeth, chipped teeth and stumps.

3. Contraindication for parapulpal retention pins:

The pin systems may only be screwed into healthy, caries-free dentine.

The specified speeds for the rotary instruments may not be exceeded. The parapulpal retention pins are single-use products and must not be reused. The rotary instruments can be reused. Parapulpal retention pins must not be used if less than 1mm of circular space is available in the dentine. (keep approx. 1 mm distance to the tooth enamel.)

4. Notes on processing:

The rotary instruments and the parapulpal retention pins are supplied non-sterile and must be cleaned, disinfected and sterilised prior to initial use. Place the rotary instruments into a cleaning and disinfection bath immediately after use. This is the only way to prevent drying in of the tissue substances. Then clean, check, disinfect and sterilise. Use only approved aldehyde-free disinfectants with proven efficacy (e.g. VAH or FDA/EPA approval or CE labelling).

Attention: Never use alcohol-based disinfectants! Additional ultrasonic cleaning (only in the bur block to avoid contact damage) is potentially possible. Observe the concentration and immersion time according to the manufacturer's instructions for the disinfectant bath. All instruments can also be cleaned in a thermal disinfectant (use bur block 4063) (other than colour-anodised mandrels).

Inspecting the instruments:

Damaged, bent, worn and non-concentric instruments are to be sorted out immediately and no longer to be used.

Then seal the instruments in foil (disposable sterilisation packaging) and sterilise in an autoclave (Class B) (134°C / 273° F; 2.1bar; 5min).

Detailed processing recommendations according to DIN EN ISO 17664

➔ see www.nti.de

5. Recommendations for use

The use of a rubber dam is recommended as a matter of principle. The dentist is responsible for selecting the appropriate retention pin. Before placing the composite filling, we recommend covering the pen with colour.

5.1. Retopin Gold, Reto-Minipin Gold and Retopin Titanium, Reto-Minipin Titanium

- ① Use the appropriate pre-drill (observe the colour marking) to drill in up to the stop and then out in a single swift movement.
Speed: 1000 rpm. Repeated drilling in and out will widen the hole too much.
- ② Use an contra-angle to screw in the pin at max. 750 rpm or with the blue handle. The pin will break automatically at the predetermined breaking point.
- ③ If several pins are required to anchor the filling, the process is repeated. In some cases the pin has to be bent for the following build-up.
The Reto-Minipin cannot be bent. Then proceed with the filling or build-up according to the manufacturer's instructions.
The mandrels can be reused. Unscrew the remaining part of the pin with pliers and replace with a new pin. Perform processing as described in Section 4.

5.2. TRI-STAR

- ① Use the appropriate pre-drill (observe the colour marking) to drill in up to the stop and then out in a single swift movement.
Speed: 1000 rpm. Repeated drilling in and out will widen the hole too much.
- ② Use an contra-angle to screw in the pin at max. 750 rpm or with the blue handle. The pin will break automatically at the predetermined breaking point.
- ③ If several pins are required to anchor the filling, the process is repeated. Then proceed with the filling or build-up according to the manufacturer's instructions. The mandrels are not intended for reuse and are disposed of.

5.3. TRI-JET

Note: the conical shank of the TRI-JET ensures that the retention pins run slightly out of round in the contra-angle.

This makes it easier to screw in the TRI-JET in poorly visible places.

- ① Use the appropriate pre-drill (observe the colour marking) to drill in up to the stop and then out in a single swift movement.
Speed: 1000 rpm. Repeated drilling in and out will widen the hole too much.
- ② Use an contra-angle to screw in the pin at max. 750 rpm or with the blue handle. The pin will break automatically at the predetermined breaking point.
- ③ Then proceed with the filling or build-up according to the manufacturer's instructions. The mandrels are not intended for reuse and are disposed of.

Storage and Keeping Conditions:

Keep dry and protected against recontamination. The instruments are to be protected against chemicals, acids, heat and extreme temperature fluctuations as a matter of principle.

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