NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
e-mail: nti@nti.de • www.nti.de

DE Produktinformation FlexiPol Polierscheiben

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← www.nti.de

GB Product information FlexiPol polishing discs

Areas of application: Separating, contouring, finishing and high-lustre polishing of all current composites.

Material composition: Stainless steel mandrel.
Polisher with silicone section for connecting to the mandrel.
Flexible plastic disc with integrated abrasives.

Requirements: The composite restoration should be prepared and smoothed prior to polishing with tungsten carbide finishers or with fine or superfine diamond instruments. The interdental region can be prepared and contoured with a fine or superfine diamond finishing strip.

Using FlexiPol polishing discs:

The colour of the discs indicates the degree of abrasiveness.

Green - single-sided, coated on top surface: initial shaping, contouring.
Blue - coated on both sides: pre-polishing, removal of scratches.
Red - coated on both sides: fine polishing.
White - coated on both sides: high-lustre polishing.

Note: FlexiPol polishing discs are single-use products.

FlexiPol polishing discs should be disinfected before initial use. Mandrels can be reused following preparation. FlexiPol polishing discs should only be used with the M4038 mandrel on patients.

FlexiPol discs are used wet (min. 50ml/minute). Polish using light, circular movements on the composite surface. The optimal motor speed is 10,000 - 12,000 rpm for approx. 10 to 20 seconds. The finer the grit size, the lower the motor speed. Excessive polishing pressure (do not exceed max. 3 N) and too high a motor speed result in overheating and damage of the composite. Avoid contact with the gingiva to safeguard against the risk of injury. Jamming or excessive bending damages the FlexiPol polishing discs.

Light polishing pressure and circular movements facilitate attaining optimum polishing results. The polishing sequence is green, blue, red, white.

Preparing, checking, disinfecting and sterilizing FlexiPol polishing discs:

FlexiPol polishing discs are single-use products. The discs should be disinfected before using for the first time. FlexiPol polishing discs cannot be sterilized. Discard polishing discs after use!

Place the mandrels in a disinfectant and cleaning solution approved for rotary instruments before initial use and after each use. The mandrels should then be checked for defects under a

magnifying glass with ten times magnification. Dispose of any defect mandrels and then sterilize undamaged mandrels. Mandrels should be sterilized in an autoclave at 134°C / 5 minutes.

CAUTION: The disinfectant solution should not contain alcohol and must be approved for use with polishers. Solutions that contain alcohol damage polishers and this may result in injury to the operator or patient.

Detailed preparation recommendations according to DIN EN ISO 17664 ← www.nti.de

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