NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • www.nti.de

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(GB) Instructions for use

FiberMaster TopHead special heads made from glass-fibre for use with conical FiberMaster glass-fibre root posts.

Contents:

1. **Product description**
2. **Indications of the FiberMaster TopHead**
3. **Contraindications**
4. **Preparation instructions**
5. **Requirements for use of the conical FiberMaster posts**
6. **Recommended uses**

6.1. Application as a post system with head

6.2. Application for temporary restorations without a post

1. Product description:

Dear customer, TopHead root post heads for conical FiberMaster glass-fibre posts should only be also used by dentists and fitted in patients. TopHeads are a system, which can be used independently of the conical FiberMaster.

Please read the product information carefully before using TopHead special heads for the conical FiberMaster root post system. FiberMaster posts are made of a material specially developed for use in the dental sector and should be used as recommended. The manufacturer does not accept any liability for damages, caused by failure to follow these instructions or use that is not compatible with the system. The user is obliged to accept responsibility for checking that the material is suitable for the intended use. This is particularly required if it involves use which is not named in the product information.

2. Indications of the FiberMaster TopHead

With coronally severely damaged or fractured teeth where the roots have no further lesions and are still fully intact periodontally. TopHeads can also be used for retaining temporary crowns. This enables endodontic treatment to be continued without impairing the aesthetics.

3. Contraindications:

TopHeads are not suitable for use with teeth with fractured or cracked roots. TopHeads are not suitable for use with teeth, which have dentine that is not sufficiently stable for anchorage of a root post.

4. Preparation instructions:

The rotary instruments should be disinfected in the usual way, sterilised in an autoclave (134°C/273°F, 5 minutes) and re-used.

The posts can be disinfected using 70% medical alcohol and any grease cleaned off. Sterilisation is not possible. Do not touch with bare fingers after cleaning.

5. Requirements for use of conical FiberMaster posts:

Endodontic treatment must be successfully completed. The root must be periodontally intact and caries-free. The post of the conical FiberMaster should contact the root wall for min. 2/3 the length of the canal. Adequate root wall thickness (min. 1.5 mm). Select as large a post diameter as possible. The occlusal or incisal surface of the tooth preparation must be prepared for a core.

A 1.5 - 2.0 mm ferrule preparation in the natural tooth structure improves the stability. Ferrule preparation to this depth is not absolutely essential due to the unique design of the FiberMaster TopHead. All instruments, posts and heads of the FiberMaster TopHead system should be secured against aspiration and swallowing. The use of a rubber dam is recommended. All instruments of the FiberMaster Systems should always be prepared before use!

6. Recommended uses:

6.1. Application as post system with head

FiberMaster posts are intended for single use and should not be re-used. Smooth the surface of the tooth using a diamond instrument. Select the size of the TopHead and root post with the aid of the X-ray image. The largest possible size should be selected. Please note, always use the matching posts and heads. Select the drilling depth for the conical root post individually, so that approx. 3 - 4 mm of filled canal remains in the apical region.

Several posts with and without TopHead can be used for teeth with multiple roots, depending on the available space. Remove the root filling. Mark the drilling depth of the drill of the conical FiberMaster using the universal depth gauge and extend the canal to the necessary depth (speed: opt. 500 - 1,000 rpm). If the gauge is too loose, simply activate it with pliers. Prepare a rest in the tooth surface for the TopHead using the two-tier hand grinder. To ensure the stability of the system, drill to the end of the diamond coating with the grinder (speed opt. 1,000 - 1,500 rpm).

Remove the conical FiberMaster from the pack using tweezers. Do not handle the posts, as the skin oils prevent a reliable bond. The posts are radiopaque, so the fit can be checked on an X-ray. If the posts have been handled, remove the contamination using 70% medical alcohol. Do not handle after cleaning. The posts cannot be sterilised, as the temperatures damage the material.

An adhesive luting technique and ETCH AND RINSE or Total Etch should be used. A rubber dam should always be placed for this technique.

Self-etching one-bottle systems are unsuitable for luting glass-fibres ①.

Also use tweezers to remove the head selected.

- Etch the area for the TopHead and surface of the prepared tooth using

- 37% phosphoric acid
- Rinse off the phosphoric acid
- Dry with paper tips.

- Condition the root canal and tooth surface with primer and bonder according to the manufacturer's instructions and light cure, if necessary.

The root post and TopHead can be silanised, if desired, though this is not absolutely necessary.

The luting composite must be applied to the root post and TopHead. Use a low viscosity dual composite for luting. High viscosity dual composites should be used with an ultrasonic unit. When luting, apply light pressure to the occlusal of the post and TopHead using a round instrument.

Use only dual composites, which also cure without light.

Adhere strictly to the adhesive manufacturer's instructions for use! Only adhesive systems should be used, which are also suitable for building up the core. Recommended systems are Luxabond and Luxacore from DMG, Panavia F2 from Kuraray. Systems, which use two different materials for luting and the head, are unsuitable for use with the TopHead system.

Use of a rubber dam is standard practice with all adhesive systems.

Caution please note:

Do not clean the canal with NaOCl (sodium hypochlorite) or H₂O₂ (hydrogen peroxide), as these chemicals can have a detrimental effect on the polymerisation of the composites. These chemicals can inhibit curing and prevent the bond in the canal. Distilled water should be used for rinsing.

Apply the adhesive to the root post and TopHead. First, insert the conical FiberMaster into the canal applying gentle pressure and then fit the TopHead. Press in the two components slowly to allow excess to escape.

Please note: It is particularly important that the head is inserted until it is in contact with the step. If required, tap the post with an instrument or use an ultrasonic round instrument to aid insertion.

Use the excess for building up the core. This creates a monoblock made from a single material, which improves the stability. Then light cure the adhesive.

Caution: The root post only transmits light apically and not sufficiently to the walls of the canal. This is why self-curing or dual-curing composite adhesives should be used.

Fabricate the core using the same material used for luting the post (please follow the manufacturer's instructions). The resulting monoblock provides greater stability than if two materials had been used. A ferrule preparation of 1,5 - 2 mm into the root surface is recommended. If the tooth is broken too far down, a deep preparation is no longer necessary.

6.2. Application for temporary restorations without a post

If there is no guarantee that endodontic treatment can be completed successfully and the endodontic treatment must be continued, the TopHead can be used for reliably retaining a temporary crown. In order to provide the temporary crown with adequate stability, the TopHead is used without the FiberMaster root post.

Select the correct size of TopHead with the aid of an X-ray. Prepare a rest in the tooth surface for the TopHead using the two-step head grinder. To ensure the stability of the system, drill to the end of the diamond coating with the grinder (speed opt. 1.000 - 1.500 rpm).

Remove the selected head using tweezers and lute adhesively as described above. Insert a paper tip, cotton pellet or a foam rubber pellet into the hole of the TopHead and root canal during luting. This keeps the access to the apex intact. Contour as required. Also leave a spacer in the TopHead before cementing the temporary crown.

Once the endodontic treatment is complete, the root post can be adhesively luted. The canal and TopHead should be cleaned and roughened with the appropriate roughening instrument to produce a reliable connection to the TopHead and canal wall.

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Detailed preparation recommendations according to DIN EN ISO 17664 www.nti.de

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© NTI-Kahla GmbH Rotary Dental Instruments