NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • http://www.nti.de

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GB Instructions for use

Titanium alloy (Ti6Al4V) post and core systems:
HeadMaster System active and passive with conical and cylindrical head
PinMaster System active

PostMaster System passive

Contents:

1. Product description
2. Indication of the post and core systems
3. Contraindication
4. Preparation instructions
5. Requirements for use of the post and core system
6. General recommended uses
- 6a. Recommended uses with HeadMaster active, PinMaster
- 6b. Recommended uses with HeadMaster passive, PostMaster

1. Product description:

For use in the dental practice by properly qualified operators (degree in dentistry). Post-and-core systems are dental precision posts and should only be used for the specific intended purpose. They are used for anchoring plastic filling materials (amalgam, composites), core build-ups in the tooth.

Please read the product information carefully before using the root post system.

The manufacturer does not accept any liability for damages caused by non-compliance with these guidelines or non-system use.

The user is obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if the application is not stated in the product information.

2. Indication of the post-and-core systems

Teeth with heavily damaged coronal sections or fractured teeth with roots with no further lesions and which are still fully intact periodontally.

3. Contraindication of the post-and-core systems

They are unsuitable for use with teeth whose roots are fractured or cracked. They are unsuitable with teeth whose dentine is no longer stable enough to reliably anchor a post-and-core system.

4. Preparation instructions:

Before using the system, disinfect it in an autoclave, rinse and sterilise: temperature 134°C, pressure 2.1 bar, time 5 min.

Detailed preparation recommendations according to DIN EN ISO 17664 www.nti.de

5. Requirements for use of the post and core systems:

Root treatment must be successfully completed. The root must be periodontally sound and caries-free. The posts should fit flush to the wall for a minimum of 2/3 the canal length.

Adequate root wall thickness (min. 1.5 mm).

Select the largest possible pin diameter.

The remaining occlusal and incisal structure must be prepared for a core build-up.

A ferrule preparation with 1.5 - 2.0 mm in the natural tooth structure improves the stability and retention of the subsequent restoration.

All instruments and root posts of the system should be secured against aspiration and swallowing. Use of a rubber dam is recommended.

All instruments and root posts of the system should always be prepared before use!

Please note, the titanium root posts are intended for single use only. Do not reuse!

6. General recommendations for use:

The post and core systems can be inserted with flowable phosphate cement. The following steps are performed the same with all post and core systems.

1. Remove the root filling.
 2. Prepare the root using a reamer.
 3. Smooth the root surface, for example using a round diamond bur.
 4. Select the size required using an X-ray.
 5. Only required with posts with a head (HeadMaster): Drill to a depth of approx. 1 mm using a facing cutter (speed: 2,000 min⁻¹) ①⑪.
 6. Determine the required drill depth using the X-ray and mark with an end-stopper.
 7. Drill the location hole to the planned depth (speed: 500-1,000 min⁻¹) ②⑥⑫⑬.
 8. Clean and disinfect the root canal.
 9. Attach the post to the insertion device.
- 6a. Procedure with:**
active HeadMaster posts with head;
PinMaster, active without head:
10. First screw the post carefully into the canal without cement ③⑦. After two screw turns, unscrew the post one turn. This prevents jamming, reducing the risk of root fracture. Screw in to the final position (with HeadMaster it is essential to insert it until the head lies on the tooth ④) and then unscrew.
 11. Apply cement to the shank of the post.

12. Fill the canal with cement.
13. Carefully insert the active post in the canal.
14. After curing, the post can be additionally secured against twisting using a parapulpal post ⑨.
15. Fabricate the build-up using plastic filling material ⑤⑩.

6b. Procedure with:

passive HeadMaster posts with head;

PostMaster, passive without head:

10. Fill the canal with cement.
11. Apply cement to the shank of the post.
12. Insert the post in the canal applying light pressure and tap in to the apex (with HeadMaster it is essential to insert it until the head lies on the tooth ⑭).
13. After curing, the post can be additionally secured against twisting using a parapulpal post ⑩.
14. Fabricate the build-up using plastic filling material ⑤⑩.

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© NTI-Kahla GmbH Rotary Dental Instruments