NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
nti@nti.de • www.nti.de

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GB Instructions for use for NTI Endo Clean instruments

Contents:

1. Product description
2. Indication of NTI Endo Clean instruments
3. Contraindication of the NTI Endo Clean instruments
4. Preparation instructions
5. Recommended uses and treatment protocol

1. Product description:

Dear customers

NTI EndoClean instruments are made from nickel titanium.

Please read through the instructions for use carefully before using NTI Endo Clean instruments.

NTI Endo Clean instruments were developed for use by qualified dentists and may only be used by them. The user is obliged to accept responsibility for checking that the material is suitable for the intended use. This is especially necessary if it involves a use not mentioned in the product information. The manufacturer does not accept liability for any damage due to the user disregarding the guidelines of these instructions or use that is not compatible with the system

2. Indication of NTI Endo Clean instruments:

NTI Endo Clean instruments were developed for endodontic use in root canal treatment.

3. Contraindication of NTI Endo Clean instruments:

NTI Endo Clean instruments should never be used for any application not listed in the methods below. In particular, the torques and speeds listed should not be exceeded.

4. Preparation instructions:

Damaged, bent or instruments that no longer run concentrically should be immediately removed and no longer used. It is recommended only to use nickel titanium instruments once. This reduces the risk of fracture.

Gates-Glidden instruments should be disinfected and sterilised in the usual way: use aldehyde-free disinfectants, VAH or FDA approved and CE Mark, possibly with an ultrasonic cleaner (only in stands), concentration according to the manufacturer's instructions. Gates-Glidden instruments can also be cleaned in a thermodesinfektor, then shrink-wrapped and sterilised in an autoclave (134°C / 273°F; 2.1bar) and re-used.

Adhere to the Robert-Koch Institute guidelines (hygiene requirements in dentistry) and the valid legal and hygiene regulations in your region!

Detailed preparation recommendation pursuant to DIN EN ISO 17664

☛ www.nti.de

5. Recommended uses:

NTI Endo Clean instruments may only be used with an endodontic contra-angle with torque control.

Avoid use at an unsuitable angle, inadequate straight access or applying excessive pressure (increased risk of fracture). Operational rpm and torque guidelines: the operational rpm and torque settings are listed in the respective preparation sequence.

Electrical measurement of the length using an Apex Finder and/or X-ray image is absolutely essential.

Every root canal treatment must be performed using a rubber dam.

5.1. Sequence 1: conventional root canal treatment

Create a glide path to the full working length manually using a K-file ISO 015 extending to the apex in the balanced force technique.

Information:

Use the balanced force technique with K-files (REF: 173.25.654.015). These should be inserted into the root canal applying light apical pressure around 90-180° in a clockwise direction. Remove hard structure at the final rotation of 120-180° in a counter clockwise direction while maintaining the pressure apically; the file should then be removed in a clockwise direction. (Color Atlases of Dental Medicine, Endodontology: Michael A. Baumann; Rudolf Beer).

After creation of the glide path, ream the canal to the first coronal curvature (5 Ncm, 600 min⁻¹) using a Gates Size 1.

Then ream the canal using a Gates Size 2 to 2mm before the first coronal curvature (5 Ncm, 600 min⁻¹).

Then ream the canal using a Gates Size 3 to 3mm before the first coronal curvature (5 Ncm, 600 min⁻¹).

Prepare to the full working length using the nickel-titanium file 20/.04 (REF: NT04-020, marking two rings and a yellow coloured ring) and 1.5 Ncm with 300 min⁻¹.

Prepare to the full working length using the nickel-titanium file 25/.04 (REF: NT04-025, marking two rings and a red coloured ring) and 1.6 Ncm with 300 min⁻¹.

With wide canals now prepare to the full working length using the nickel-tita-

nium file 35/.04 (REF: NT04-035, marking two rings and a green coloured ring) and 1.8 Ncm with 300 min⁻¹.

Then fill the root canal to the apical constriction and finally a radiological check.

5.2. Séquence 2: minimally invasive root canal treatment

Create a guide path to the full working length manually using a K-file ISO 015 (REF: 173.25.654.015) extending to the apex in the balanced force technique.

Prepare the access extension to the first coronal curvature using the nickel-titanium file 25/.06 (REF: NT06-025, marking three rings and a red coloured ring) and 2.5 Ncm with 600 min⁻¹.

Prepare to the full working length using the nickel-titanium file 25/.04 (REF: NT04-025, marking two rings and a red coloured ring) and 1.6 Ncm with 300 min⁻¹.

With wide canals now prepare to the full working length using the nickel-titanium file 35/.04 (REF: NT04-035, marking two rings and a green coloured ring) and 1.8 Ncm with 300 min⁻¹.

Then fill the root canal to the apical constriction and finally a radiological check.

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