NTI-Kahla GmbH • Rotary Dental Instruments
 Im Camisch 3 • D-07768 Kahla/Germany
 Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
 e-mail: nti@nti.de • http://www.nti.de

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GB Preparation of Bur Stands in Dental Practices

- 1) Place the bur stand with the instruments in a disinfectant approved for rotating instruments (can possibly also be used for polishers, depending on contents). Never use disinfectants containing alcohol. It is possible to carry out a thermal disinfection.
- 2) Remove the burs after disinfecting, rinse each individual stand and blow dry. Ensure that every hole has been blown dry separately.

3. Rearrange the instruments afterwards.
4. Because of its low height, the bur stand can be shrink-wrapped in an autoclave pouch.
5. It is then sterilised in an autoclave at 134 °C (273 °F) 5 minutes or 121 °C (250 °F) 20 minutes.
6. Finally store the shrink-wrapped stands along with the instruments at room temperature in dust-free conditions and protected from direct sunlight (to protect the polishers).

Manual cleaning and disinfection			Preparation for sterilisation	Sterilisation	Sterilisation approval und storage
Process	Substance	Exposure time			
Chemical cleaning and disinfection ¹⁾	Aldehyde-free disinfectant, VAH or FDA approved and CE mark, may be used with ultrasonic cleaner, concentrations according to the manufacturer's instructions.	according to manufacturer's instructions	Remove residual disinfectant with water and dry ²⁾ Packaging complies with EN ISO 11607/DIN EN 868-2ff.	Steam sterilization 121°C / 20 min. 134°C / 5 min. Steam sterilizer should comply with DIN EN 13060 or DIN EN 285.	Approval or rejection in the case of incorrect procedure to be carried out by specialist Storage: In a dry place and protected from contaminants

Detailed recommendations for preparation in accordance with DIN EN ISO 17664 www.nti.de

¹⁾Use only disinfectants and cleaning agents which are suitable for aluminium. Do not use alkaline bur baths. Avoid the use of disinfectants containing alcohol! ²⁾The use of inappropriate disinfectants and/or leaving the instruments in the disinfectant too long can result in corrosion, discolouration or increased chance of damage occurring.

Observe the valid legal and hygiene regulations in your country!

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